

CHEF GEL WITH EMBEDDED YEAST DNA

- 1) Inoculate a single colony into 50 ml YPD. Grow in shaker to an OD₆₀₀ of > 1.0 at 30° C.
- 2) Turn on 50° C water bath. When the desired OD is reached, centrifuge the cells at 5,000 x g, 10 minutes at 4° C. Decant the supernatant and resuspend the pellet in 10 ml cold 50 mM EDTA pH 8.0.
- 3) Determine the cell concentration by direct cell counts. 1:100 dilution usually works. You want 6×10^8 cells/ml of agarose plugs (3×10^8 cells per 0.5 ml plugs).
- 4) Microwave the 2% CleanCut agarose solution to melt, equilibrate to 50° C.
- 5) Remove appropriate number of cells for amount of plugs desired. Resuspend the cells in 500 μ l (250 μ l) of cell suspension buffer and bring suspension to 50° C. (Numbers in parentheses refers to the amount for 0.5 μ l of plugs).
- 6) Just prior to mixing cells with agarose add 50 μ l (25 μ l) of 20 mg/ml Lyticase and proceed immediately to next step.
- 7) Combine warmed cell suspension with 500 μ l (250 μ l) of 2% CleanCut agarose and mix gently but thoroughly. Transfer the mixture to plug molds. Allow the agarose to solidify- 15 minutes at 4° C.
- 8) In a Falcon tube add 5 ml (2.5 ml) of lyticase buffer per strain/plugs. Push the solidified plugs, using the provided tool, into the Falcon tube with lyticase buffer. Incubate the plugs 1 hour at 37° C without agitation.
- 9) Remove the lyticase buffer and rinse the plugs with ddH₂O. Remove ddH₂O completely and add 5 ml (2.5 ml) of Proteinase K reaction buffer. Incubate the plugs overnight at 50° C without agitation.
- 10) Wash the plugs 4x in 30-50 ml of wash buffer, 1 hour per wash. Wash at room temperature with gentle agitation. Wash with 1 mM PMSF during the second or third wash if the plugs are to be used in subsequent enzymatic reactions.
- 11) Store the plugs at 4° C. The plugs are stable for 3 months to 1 year. For long term storage, store in wash buffer.

Cell Suspension Buffer

10 mM Tris, pH 7.2

20 mM NaCl

50 mM EDTA

Proteinase K Reaction Buffer

100 mM EDTA, pH 8.0

0.2% Sodium deoxycholate

1% Sodium lauryl sarcosine

1 mg/ml Proteinase K

Lyticase Buffer

10 mM Tris, pH 7.2

50 mM EDTA

1 mg/ml lyticase

1X Wash Buffer

20 mM Tris, pH 8.0

50 mM EDTA