

Denaturing purification of Insoluble Proteins

<u>Buffer A</u>	<u>500 ml</u>	<u>1 L</u>
6 M Guanidine hydrochloride	286.59 g	573.18 g
0.1 M NaH ₂ PO ₄	6.9 g	13.8 g
0.01 M Tris	5 ml	10 ml (of 1 M Tris HCl pH 8.0)
adjust pH to 8.0 with NaOH		

<u>Buffer B</u>		
8 M Urea	240.24 g	480.48 g
0.1 M NaH ₂ PO ₄	6.9 g	13.8 g
0.01 M Tris	5 ml	10 ml (of 1 M Tris HCl pH 8.0)
adjust pH to 8.0 with NaOH		

<u>Buffer C</u>		
8 M Urea	240.24 g	480.48 g
0.1 M NaH ₂ PO ₄	6.9 g	13.8 g
0.01 M Tris	5 ml	10 ml (of 1 M Tris HCl pH 8.0)
adjust pH to <u>6.3</u> with HCl		

<u>Buffer E</u>		
8 M Urea	240.24 g	
0.1 M NaH ₂ PO ₄	6.9 g	
0.01 M Tris	5 ml (of 1 M Tris HCl pH 8.0)	
adjust pH to <u>4.5</u> with HCl		

1. Resuspend the cells in Buffer A: 5 ml per gram weight.
2. Stir cells for 2 hours at room temperature.
3. Centrifuge lysate at 10,000 rpm for 30 minutes at 4°C, collect supernatant.
4. Take 1 ml of a 50% slurry of Ni-NTA resin, wash 3 times in Buffer A (50 ml), equilibrate in Buffer A (v/v).
5. Add 50% Ni-NTA resin to supernatant from step 3 and stir at room temperature for 1-2 hours.
6. Centrifuge at 1600 rpm for 5 minutes. Discard supernatant.
7. Wash beads 3 times with 10-50 ml of Buffer A (centrifuge at 1600 rpm for 5 minutes).
8. Wash beads 3 times with 10-5- ml of Buffer B.
9. Wash beads 3 times with 10-50 ml of Buffer C.
10. Elute protein with 500 μ l of Buffer E (7 times).
11. Centrifuge the eluates for 20 seconds at 10,000 rpm to remove rest of beads.
12. Analyze samples on SDS-PAGE: 20 μ l sample + 10 μ l 3x Lämmli buffer.

Protein gel electrophoresis (Lämmli)

- solutions:
- Tris gel buffer: 375 mM Tris pH 8.8
0.1% SDS
 - stacking gel (4%): 26.67 ml acrylamide-stock (30%)
0.2 g SDS (0.1%)
3 g Tris base
adjust pH to 6.8
add H₂O to 200 ml
store at 4°C
 - 10x running buffer: 250 mM Tris 30.3 g/l Tris base
1920 mM glycine 144.2 g/l
1% SDS 10 g/l
pH should be ~8.3 (do not adjust pH!!!)
 - 3x Lämmli: 187.5 mM Tris pH 6.8
6% SDS
30% glycerol
15% β -mercaptoethanol
0.003% bromphenolblue
make 1 ml aliquots
heat 5 min at 100°C in heating block
- sample preparation:
- x μ l protein extract
 - 5 μ l 3x Lämmli
 - H₂O to 15 μ l
 - denature 5 min at 95-100°C (heating block)
 - samples can be stored at -20°C (heat again before loading)
- gel preparation:
- for Hoefer Minigels (Mighty Small II SE250)
 - clean 1 glass plate, 1 white plate and 2 spacers (0.75 mm) per gel
 - assemble for pouring (white plate in the back, do not tighten screws too much, make sure spacers and plates are level at the bottom)
 - prepare 10% APS (ammonium persulfate) in H₂O
 - make a mark on glass plate to indicate the volume of the gel (~1 cm below the comb)
 - acrylamide-solution: e.g. for two 10% gels:
 - 3.3 ml acrylamide-stock (30%)
 - 6.7 ml Tris gel buffer
 - 200 μ l 10% APS
 - 6 μ l TEMED
 - mix well and pour gel up to mark
 - add 200 μ l n-butanol on top
 - polymerize for ~1 h

- remove gel from pouring unit, remove n-butanol with H₂O
- remove drops between plates using a whatman paper
- stacking gel: for 2 gels:

4 ml stacking gel (4%)
 200 μ l 10% APS
 4 μ l TEMED

- mix well and fill gel up to the very top
- immediately add combs (stacking gel polymerize very fast)
- polymerize for 15 min

-gel run:

- prepare 250 ml 1x running buffer
- carefully remove combs and load samples
- add protein standard
- run gels for 50-60 min at 25-30 mA per gel (use water cooling and mix buffer) until bromphenolblue reaches bottom of gel and runs out
- the protein gels can be stained directly or used in a western blot procedure

Staining of protein gels with Coomassie

- solutions:
 - Coomassie-solution 1: 0.025% Coomassie Brilliant Blue R250
25% isopropanol
10% HAc
 - Coomassie-solution 2: 0.0025% Coomassie Brilliant Blue R250
10% isopropanol
10% HAc
 - destaining solution: 10% isopropanol
10% HAc

- shake protein gel 20-30 min in Coomassie-solution 1
- partially destain in Coomassie-solution 2 (gels can be left in this solution almost indefinitely)
- destain 2-3 h in destaining solution (add a rolled kim wipe to absorb the coomassie)
- fix for ~1 h in 10% HAc
- sandwich gel between wet cellophane foil, put sandwich on a whatman filter and cover with saran wrap
- dry for 1 h on gel drier at 80°C
- make sure gels do not curl