UC Davis γ-ray source

Contact:

Don Ballard Tel. 3-1245 pager 916-523-0037 dwballard@ucdavis.edu Principal Lab Mechanic Institute of Toxicology & Environmental Health 120 ITEH

or

Shane McCarty 2-1342

Before experiment contact either Don or Shane at least 24 hr before to sign up. They run the machine for us.

Use is \$50 per hr.

After 10 am is good, pm times are best.

There is a lab bench available for post-irradiation workups like plating. There is no gas hookup.

Source:

840 R/min equivalent to 8 Gy/min ^{137}Cs source exposure chamber revolves at 6 rpm exposure chamber is about $8 \times 10 \times 12$ inches

Procedure:

- Suspend cells appropriately in water with a final volume of 500 μ l in eppendorf tube.
- Place eppendorf tubes in ice.
- Irradiate with the appropriate dose at center of plate to ensure even dose.
- Dilute in water and plate.