

***E. coli* competent cells and transformation (Inoue et al. 1990 Gene 96, 23-28)**

How to prepare TB buffer?

1. Prepare 10 mM Hepes (or Pipes)/15 mM CaCl₂/250 mM KCl
2. Adjust pH to 6.7 with KOH.
3. Dissolve MnCl₂ to have final concentration of 55 mM in TB buffer.
4. Sterilize solution by filtration through 45 µm filter and store at +4°C.

Note: All salts should be added as solids.

How to prepare competent *E. coli* (e.g. DH5[]) cells?

1. Inoculate 250 ml of SOB medium in a 2 liter flask and grow to OD₆₀₀ =0.6 at 18°C (or at room temperature) with vigorous shaking (200-250 rpm).
2. Put flask on ice for 10 min.
3. Spin down cells at 2,500 x g (*i.e.* 3K rpm Beckman J-6B centrifuge) for 10 min at +4°C.
4. Resuspend pellet in 80 ml of ice-cold TB and incubate for 10 min, then spin down cells.
5. Resuspend pellet gently in 20 ml of TB and add DMSO to final concentration of 7% with gentle swirling.
6. Incubate on ice for 10 min.
7. Dispense cell suspension in 0.5 ml aliquots and immediately flash-freeze by immersion in liquid N₂.
8. Store frozen cells at -70/80°C.

How to perform a transformation?

1. Thaw competent cells on ice.
2. Dispense 100 µl into 1.5 ml microcentrifuge tubes.
3. Add 0.5-5 µl of plasmid DNA in each tube and incubate on ice for 30 minutes.
4. Heat shock cells without agitation at 42°C for 90 seconds and transfer to ice bath.
5. Add 1 ml of SOC medium (or LB+) and shake vigorously at 37°C for 1 hour (tape tubes to rotating wheel).
6. Plate cells onto LB+Amp plates (50 µg/ml Ampicillin) and incubate them at 37°C over night.

SOB

0.5 % yeast extract
2 % tryptone
10 mM NaCl
2.5 mM KCl
10 mM MgCl₂
10 mM MgSO₄

SOC

0.5 % yeast extract
2 % tryptone
10 mM NaCl
2.5 mM KCl
10 mM MgCl₂
10 mM MgSO₄
20 mM glucose

LB+

5 ml LB
50 µl 40% glucose
100 µl 1M MgCl₂

Revised: 2/06