Restriction analysis of plasmids

- Solutions: - 1x TBE: make from 10xTBE

10xTBE: 0.89M Tris base 109 g

 $0.89M H_3BO_3$ 55 g 0.025M EDTA 9.3 g

add H₂O to 11

- 10x loading buffer: 20% Ficol 1400 10 g

 0.1 M EDTA
 10 ml 0.5 M

 1% SDS
 5 ml 10%

 0.25% bromphenolblue
 0.125 g

 0.25% xylencyanol
 0.125 g

 add H₂O to 50 ml

- use ~200 ng plasmid DNA

- 10 μl sample: - 1 μl 10x buffer

- x µl DNA

- x µl restriction enzyme (1 U cuts 1 µg DNA in 1 h at 37°C)

-(0.1 µl BSA 10 mg/ml if needed)

- add H_2O to $10 \mu l$

- incubate 1 h at 37°C
- add 1 µl 10x loading buffer
- 0.8% agarose gel: dissolve 0.8 g agarose in 100 ml TBE-buffer (heat in microwave)

- add 5 μ l 10 mg/ml ethidium bromide (= 0.5 μ g/ml)

(wear gloves, EtBr is carcinogenic!)

- pour into gel box and let polymerize (~45 min)

- load samples on gel and run at 5-10 V/cm length of gel
- don't forget to add 5 µl of DNA ladder as standard
- when blue markers are sufficiently separated, take a photograph on UV table