

SDS-PAGE Protocol

SDS-PAGE Solutions

40% Acrylamide (37.5:1)

Acrylamide	116.8 g
N,N'-Methylene bisacrylamide	3.2 g
DDI H ₂ O	to 300 ml

Filter and store in a dark bottle at 4°C. (We buy this premade)

30% Ammonium Persulfate

Ammonium Persulfate	1.5 g
DDI H ₂ O	5 ml

Store at 4°C. Replace every month.

RG Buffer--1.5 M Tris•Cl, pH 8.8

DDI H ₂ O	300 ml
Tris-free base	90.75 g
Conc. HCl	8 ml

Adjust to pH 8.8 with conc. HCl, and bring final volume to 500 ml with DDI H₂O.

SG Buffer--1.0 M Tris•Cl, pH 6.8

DDI H ₂ O	300 ml
Tris-free base	60.54 g
Conc. HCl	36 ml

Adjust to pH 6.8 with conc. HCl, and bring final volume to 500 ml with DDI H₂O.

4x SDS-PAGE Sample Buffer

125 mM Tris•HCl, pH 6.8	1 M	5 ml
20% Glycerol		8 ml
4% SDS	20%	8 ml
10% β-Mercaptoethanol		4 ml
0.5 mg/ml Bromophenol Blue		20 mg
<u>DDI H₂O</u>		<u>15 ml</u>
Total		40 ml

10x SDS-PAGE Running Buffer

30.3 g	Tris base
144.0 g	Glycine
10.0 g	SDS

Dissolve and bring total volume to 1,000 ml with deionized water. Do not adjust pH with acid or base (pH is normally 8.3 as prepared).

Coomassie Stain Solution

Ethanol	150 ml
Glacial Acetic Acid	50 ml
DDI H ₂ O	300 ml
Coomasie Brilliant Blue-R-250	1 g

Dissolve Coomasie Brilliant Blue-R-250 in EtOH first.

Destain Solution

Ethanol	1200 ml
Glacial Acetic Acid	400 ml
DDI H ₂ O	2.4 l

Mini-Protean SDS-PAGE Protocol

Casting the Gel

- 1] Assemble glass plates and spacers in gel casting apparatus—see BioRad instruction manual.
- 2] Mix the components for the resolving gel as described in the Mini-Protean II protocol.
- 3] Pour the resolving gel mixture into the gel plates to a level 2 cm below the top of the shorter plate.
- 4] Pace a layer of DDI H₂O over the top of the resolving gel to prevent meniscus formation in the resolving gel.
- 5] Allow resolving gel to stand 30 min at room temperature.
- 6] Drain the DDI H₂O from top of the resolving gel. Rinse with DDI H₂O, drain, and wick any remaining DDI H₂O away with a Kimwipe.
- 7] Mix components for stacking gel.
- 8] Pour stacking gel solution into gel plates (on top of running gel), so that gel plates are filled. Insert comb to the top of the spacers.
- 9] Allow gel to stand for at least 1 hr at room temperature, or overnight at 4°C (wrapped in saran wrap).

Preparing Samples

Note: 10 well combs will hold up to 30 μ l of prepared sample. 15 well combs will hold up to 20 μ l of prepared sample.

Cell Samples

- 1] Harvest 100 μ l of cells at O.D. > 0.6. Decant the supernatant media.
- 2] Resuspend cells in 20 μ l of 2x sample buffer.
- 3] Incubate tubes in boiling water for 5 min.
- 4] Centrifuge at 12,000 x g for 30 s.

Solution Samples

- 1] Place a volume of protein solution (or 1 μ l of standard) into a μ fuge tube, such that there

is 5-10 μg of protein in the solution.

- 2] And an equal volume of 2x sample buffer (or 10 μl for standards).
- 3] Incubate tubes in boiling water for 5 min.

- 4] Centrifuge at 12,000 x g for 30 s.

Running the Gel

- 1] Remove comb and assemble cast gel into Mini-Protean II apparatus.
- 2] Add freshly prepared 1x running buffer (300 ml) to both chambers of the apparatus.
- 3] Load the prepared samples into the wells of the gel.
- 4] Run the gel at 100 V until the dye front migrates into the running gel (~15 min), and increase to 200 V until the dye front reaches the bottom of the gel (~45 min).

Staining & Destaining the Gel

- 1] Remove the run gel from the apparatus and remove the spacers and glass plates. Place the gel into a small tray. *Note:* Never use a metal spatula to separate the glass plates.
- 2] Add ~20 ml staining solution and stain for > 30 min with gentle shaking.
- 3] Pour off and save the stain.
- 4] Add ~5 ml destain solution and destain for ~1 min with gentle shaking.
- 5] Pour off and discard the destain solution. Add ~30 ml of destain solution.
- 6] Destain with gentle shaking until the gel is visibly destained (> 2 hr).
- 5] Pour off and discard the destain solution.
- 6] Rinse with DDI H₂O. Add ~30 ml DDI H₂O and rinse for 5 min with gentle shaking.
- 7] Dry the gel on the gel dryer at 60°C for 1 hr with a sheet of Whatman filter paper below the gel and a piece of Seran wrap over the gel.

% Acrylamide in running gel	Separation size range (kDa)
Single percentage:	
5%	100–250
7.5%	40–200
10%	30–150
12%	20–120
15%	10–100
18%	6–50
Gradient:	
4–15%	20–250
4–20%	10–200
10–20%	10–100
8–16%	6–70

From the BioRad Readygel Manual.

Adapted from Laemmli, U.K. (1970) Nature 227, 680-685. (see page 681).

